

KERRISON PUNCH; 40° UP CUT

5mm; with ejector; thin footplate

18.0cm (7.2"); 10mm open; CC-coated



SKU: 240-071-05ATCC

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---|
| UNIT | pcs. |
| TYPE | Laminectomy Punches |
| AUTHORS | CLEANWAVE , KERRISON |
| SPECIFICATIONS | 40° up cut , with ejector |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---|
| UNIT | pcs. |
| TYPE | Laminectomy Punches |
| AUTHORS | CLEANWAVE , KERRISON |
| SPECIFICATIONS | 40° up cut , with ejector |