

KERRISON PUNCH; 40° UP CUT

1mm; thin footplate

18.0cm (7.2"); 10mm open; CC-coated




SKU: 240-071-01TCC

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
TYPE [Laminectomy Punches](#)
AUTHORS [CLEANWAVE](#), [KERRISON](#)
SPECIFICATIONS [40° up cut](#), [without ejector](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
TYPE [Laminectomy Punches](#)
AUTHORS [CLEANWAVE](#), [KERRISON](#)
SPECIFICATIONS [40° up cut](#), [without ejector](#)