

KERRISON PUNCH; 40° UP CUT

1mm; without ejector; thin footplate

20.0cm (8"); 10mm open; CC-coated




SKU: 240-072-01TCC

Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
AUTHORS	CLEANWAVE , KERRISON
SPECIFICATIONS	40° up cut , without ejector



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
AUTHORS	CLEANWAVE , KERRISON
SPECIFICATIONS	40° up cut , without ejector