

KERRISON PUNCH; 40° UP CUT

2mm; with ejector; thin footplate

20.0cm (8"); 10mm open




SKU: 240-072-02AT

Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
MATERIAL	Ceramic(CC)
AUTHORS	CLEANWAVE , KERRISON
Medicon	57.09.70
SPECIFICATIONS	40° up cut , with ejector



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
MATERIAL	Ceramic(CC)
AUTHORS	CLEANWAVE , KERRISON
Medicon	57.09.70

SPECIFICATIONS

[40° up cut, with ejector](#)