

KERRISON PUNCH; 40° UP CUT

3mm; without ejector

20.0cm (8"); 10mm open




SKU: 240-072-03

Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|--|
| UNIT | pcs. |
| TYPE | Laminectomy Punches |
| Aesculap | FF738 |
| AUTHORS | KERRISON |
| Medicon | 57.08.83 |
| SPECIFICATIONS | 40° up cut , without ejector |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------|-------------------------------------|
| UNIT | pcs. |
| TYPE | Laminectomy Punches |
| Aesculap | FF738 |
| AUTHORS | KERRISON |
| Medicon | 57.08.83 |

SPECIFICATIONS

[40° up cut, without ejector](#)