

KERRISON PUNCH; 40° UP CUT

3mm; with ejector; thin footplate

20.0cm (8"); 10mm open



SKU: 240-072-03AT

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
MATERIAL	Ceramic(CC)
AUTHORS	CLEANWAVE , KERRISON
Medicon	57.09.71
SPECIFICATIONS	40° up cut , with ejector



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
MATERIAL	Ceramic(CC)
AUTHORS	CLEANWAVE , KERRISON
Medicon	57.09.71

SPECIFICATIONS

[40° up cut, with ejector](#)