

KERRISON PUNCH; 40° UP CUT

curved; 4mm; 10mm open;

without ejector; 20cm (8")




SKU: 240-072-04-GEB

Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
MATERIAL	Ceramic(CC)
AUTHORS	KERRISON
SPECIFICATIONS	40° up cut , without ejector



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
MATERIAL	Ceramic(CC)
AUTHORS	KERRISON
SPECIFICATIONS	40° up cut , without ejector