

KIRSCHNER WIRE DRILLS

trocar-trocar; 0.8 / 140mm

box with 10 pieces



SKU: 212-525-08

Categories: [Open Surgery](#), [Orthopaedics](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	Box
TYPE	Kirschner wires
AUTHORS	KIRSCHNER
SPECIFICATIONS	trocar-trocar

ZUSÄTZLICHE INFORMATIONEN

UNIT	Box
TYPE	Kirschner wires
AUTHORS	KIRSCHNER
SPECIFICATIONS	trocar-trocar