

## KIRSCHNER WIRE DRILLS

trocar-trocar; 1.0 / 160mm

box with 10 pieces



**SKU:** 212-526-10

**Categories:** [Open Surgery](#), [Orthopaedics](#)

### ZUSÄTZLICHE INFORMATIONEN

|                       |                                 |
|-----------------------|---------------------------------|
| <b>UNIT</b>           | <a href="#">Box</a>             |
| <b>TYPE</b>           | <a href="#">Kirschner wires</a> |
| <b>AUTHORS</b>        | <a href="#">KIRSCHNER</a>       |
| <b>Medicon</b>        | 76.10.80                        |
| <b>SPECIFICATIONS</b> | <a href="#">trocar-trocar</a>   |

### ZUSÄTZLICHE INFORMATIONEN

|                       |                                 |
|-----------------------|---------------------------------|
| <b>UNIT</b>           | <a href="#">Box</a>             |
| <b>TYPE</b>           | <a href="#">Kirschner wires</a> |
| <b>AUTHORS</b>        | <a href="#">KIRSCHNER</a>       |
| <b>Medicon</b>        | 76.10.80                        |
| <b>SPECIFICATIONS</b> | <a href="#">trocar-trocar</a>   |