

KIRSCHNER WIRE DRILLS

trocar-trocarr; 2.2 / 160mm

box with 10 pieces



SKU: 212-526-22

Categories: [Open Surgery](#), [Orthopaedics](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---------------------------------|
| UNIT | Box |
| TYPE | Kirschner wires |
| AUTHORS | KIRSCHNER |
| Medicon | 76.10.92 |
| SPECIFICATIONS | trocar-trocar |

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---------------------------------|
| UNIT | Box |
| TYPE | Kirschner wires |
| AUTHORS | KIRSCHNER |
| Medicon | 76.10.92 |
| SPECIFICATIONS | trocar-trocar |