

## KIRSCHNER WIRE DRILLS

trocar-trocar; 1.6 / 310mm

box with 10 pieces



**SKU:** 212-528-16

**Categories:** [Open Surgery](#), [Orthopaedics](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">Box</a>
<b>TYPE</b>	<a href="#">Kirschner wires</a>
<b>Aesculap</b>	LX957
<b>AUTHORS</b>	<a href="#">KIRSCHNER</a>
<b>Medicon</b>	76.13.16
<b>SPECIFICATIONS</b>	<a href="#">trocar-trocar</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">Box</a>
<b>TYPE</b>	<a href="#">Kirschner wires</a>
<b>Aesculap</b>	LX957
<b>AUTHORS</b>	<a href="#">KIRSCHNER</a>
<b>Medicon</b>	76.13.16

## SPECIFICATIONS

[trocar-trocar](#)