

## KIRSCHNER WIRE DRILLS

trocar-trocar; 1.6 / 310mm

box with 10 pieces



**SKU:** 212-528-16

**Categories:** [Open Surgery](#), [Orthopaedics](#)

### ZUSÄTZLICHE INFORMATIONEN

|                       |                                 |
|-----------------------|---------------------------------|
| <b>UNIT</b>           | <a href="#">Box</a>             |
| <b>TYPE</b>           | <a href="#">Kirschner wires</a> |
| <b>Aesculap</b>       | LX957                           |
| <b>AUTHORS</b>        | <a href="#">KIRSCHNER</a>       |
| <b>Medicon</b>        | 76.13.16                        |
| <b>SPECIFICATIONS</b> | <a href="#">trocar-trocar</a>   |

### ZUSÄTZLICHE INFORMATIONEN

|                 |                                 |
|-----------------|---------------------------------|
| <b>UNIT</b>     | <a href="#">Box</a>             |
| <b>TYPE</b>     | <a href="#">Kirschner wires</a> |
| <b>Aesculap</b> | LX957                           |
| <b>AUTHORS</b>  | <a href="#">KIRSCHNER</a>       |
| <b>Medicon</b>  | 76.13.16                        |

## SPECIFICATIONS

[trocar-trocar](#)