

KIRSCHNER WIRE DRILLS

trocar-trocar; 1.8 / 310mm

box with 10 pieces



SKU: 212-528-18

Categories: [Open Surgery](#), [Orthopaedics](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	Box
TYPE	Kirschner wires
AUTHORS	KIRSCHNER
Medicon	76.13.18
SPECIFICATIONS	trocar-trocar

ZUSÄTZLICHE INFORMATIONEN

UNIT	Box
TYPE	Kirschner wires
AUTHORS	KIRSCHNER
Medicon	76.13.18
SPECIFICATIONS	trocar-trocar