

KIRSCHNER WIRE DRILLS

trocar-round; 1.0 / 160mm

box with 10 pieces



SKU: 212-506-10

Categories: [Open Surgery](#), [Orthopaedics](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	Box
TYPE	Kirschner wires
Aesculap	LX130
AUTHORS	KIRSCHNER
Medicon	76.10.60
SPECIFICATIONS	trocar-round

ZUSÄTZLICHE INFORMATIONEN

UNIT	Box
TYPE	Kirschner wires
Aesculap	LX130
AUTHORS	KIRSCHNER
Medicon	76.10.60

SPECIFICATIONS

[trocar-round](#)