

KIRSCHNER WIRE DRILLS

trocar-round; 1.4 / 160mm

box with 10 pieces



SKU: 212-506-14

Categories: [Open Surgery](#), [Orthopaedics](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	Box
TYPE	Kirschner wires
Aesculap	LX134
AUTHORS	KIRSCHNER
Medicon	76.10.64
SPECIFICATIONS	trocar-round

ZUSÄTZLICHE INFORMATIONEN

UNIT	Box
TYPE	Kirschner wires
Aesculap	LX134
AUTHORS	KIRSCHNER
Medicon	76.10.64

SPECIFICATIONS

[trocar-round](#)