

KIRSCHNER WIRE DRILLS

trocar-round; 2.0 / 160mm

box with 10 pieces



SKU: 212-506-20

Categories: [Open Surgery](#), [Orthopaedics](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	Box
TYPE	Kirschner wires
Aesculap	LX140
KLS Martin	22-631-20
AUTHORS	KIRSCHNER
Medicon	76.10.70
SPECIFICATIONS	trocar-round

ZUSÄTZLICHE INFORMATIONEN

UNIT	Box
TYPE	Kirschner wires
Aesculap	LX140
KLS Martin	22-631-20

AUTHORS

[KIRSCHNER](#)

Medicon

76.10.70

SPECIFICATIONS

[trocar-round](#)