

KIRSCHNER WIRE DRILLS

trocar-round; 1.0 / 310mm

box with 10 pieces



SKU: 212-508-10

Categories: [Open Surgery](#), [Orthopaedics](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	Box
TYPE	Kirschner wires
Aesculap	LX030
KLS Martin	22-633-10
AUTHORS	KIRSCHNER
Medicon	76.12.60
SPECIFICATIONS	trocar-round

ZUSÄTZLICHE INFORMATIONEN

UNIT	Box
TYPE	Kirschner wires
Aesculap	LX030
KLS Martin	22-633-10

AUTHORS

[KIRSCHNER](#)

Medicon

76.12.60

SPECIFICATIONS

[trocar-round](#)