

KIRSCHNER WIRE DRILLS

trocar-round; 2.2 / 310mm

box with 10 pieces



SKU: 212-508-22

Categories: [Open Surgery](#), [Orthopaedics](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	Box
TYPE	Kirschner wires
Aesculap	LX042
KLS Martin	22-633-22
AUTHORS	KIRSCHNER
Medicon	76.12.72
SPECIFICATIONS	trocar-round

ZUSÄTZLICHE INFORMATIONEN

UNIT	Box
TYPE	Kirschner wires
Aesculap	LX042
KLS Martin	22-633-22

AUTHORS

[KIRSCHNER](#)

Medicon

76.12.72

SPECIFICATIONS

[trocar-round](#)