

KIRSCHNER WIRE DRILLS

trocar-flat; 1.5 / 160mm

box with 10 pieces



SKU: 212-516-15

Categories: [Open Surgery](#), [Orthopaedics](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	Box
TYPE	Kirschner wires
Aesculap	LX115
KLS Martin	22-630-15
AUTHORS	KIRSCHNER
Medicon	76.10.15
SPECIFICATIONS	trocar-flat

ZUSÄTZLICHE INFORMATIONEN

UNIT	Box
TYPE	Kirschner wires
Aesculap	LX115
KLS Martin	22-630-15

AUTHORS

[KIRSCHNER](#)

Medicon

76.10.15

SPECIFICATIONS

[trocar-flat](#)