

KIRSCHNER WIRE DRILLS

trocar-flat; 2.2 / 160mm

box with 10 pieces



SKU: 212-516-22

Categories: [Open Surgery](#), [Orthopaedics](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	Box
TYPE	Kirschner wires
Aesculap	LX122
KLS Martin	22-630-22
AUTHORS	KIRSCHNER
Medicon	76.10.22
SPECIFICATIONS	trocar-flat

ZUSÄTZLICHE INFORMATIONEN

UNIT	Box
TYPE	Kirschner wires
Aesculap	LX122
KLS Martin	22-630-22

AUTHORS

[KIRSCHNER](#)

Medicon

76.10.22

SPECIFICATIONS

[trocar-flat](#)