

KIRSCHNER WIRE DRILLS

trocar-flat; 2.5 / 160mm

box with 10 pieces



SKU: 212-516-25

Categories: [Open Surgery](#), [Orthopaedics](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	Box
TYPE	Kirschner wires
Aesculap	LX125
AUTHORS	KIRSCHNER
Medicon	76.10.25
SPECIFICATIONS	trocar-flat

ZUSÄTZLICHE INFORMATIONEN

UNIT	Box
TYPE	Kirschner wires
Aesculap	LX125
AUTHORS	KIRSCHNER
Medicon	76.10.25

SPECIFICATIONS

[trocar-flat](#)