

## KIRSCHNER WIRE DRILLS

trocar-flat; 1.0 / 310mm

box with 10 pieces



**SKU:** 212-518-10

**Categories:** [Open Surgery](#), [Orthopaedics](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">Box</a>
<b>TYPE</b>	<a href="#">Kirschner wires</a>
<b>Aesculap</b>	LX010
<b>KLS Martin</b>	22-632-10
<b>AUTHORS</b>	<a href="#">KIRSCHNER</a>
<b>Medicon</b>	76.15.10
<b>SPECIFICATIONS</b>	<a href="#">trocar-flat</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">Box</a>
<b>TYPE</b>	<a href="#">Kirschner wires</a>
<b>Aesculap</b>	LX010
<b>KLS Martin</b>	22-632-10

**AUTHORS**

[KIRSCHNER](#)

**Medicon**

76.15.10

**SPECIFICATIONS**

[trocar-flat](#)