

KIRSCHNER WIRE DRILLS

trocar-flat; 1.0 / 310mm

box with 10 pieces



SKU: 212-518-10

Categories: [Open Surgery](#), [Orthopaedics](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---------------------------------|
| UNIT | Box |
| TYPE | Kirschner wires |
| Aesculap | LX010 |
| KLS Martin | 22-632-10 |
| AUTHORS | KIRSCHNER |
| Medicon | 76.15.10 |
| SPECIFICATIONS | trocar-flat |

ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|---------------------------------|
| UNIT | Box |
| TYPE | Kirschner wires |
| Aesculap | LX010 |
| KLS Martin | 22-632-10 |

AUTHORS

[KIRSCHNER](#)

Medicon

76.15.10

SPECIFICATIONS

[trocar-flat](#)