

KIRSCHNER WIRE DRILLS

trocar-flat; 1.4 / 310mm

box with 10 pieces



SKU: 212-518-14

Categories: [Open Surgery](#), [Orthopaedics](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	Box
TYPE	Kirschner wires
Aesculap	LX014
KLS Martin	22-632-14
AUTHORS	KIRSCHNER
Medicon	76.15.14
SPECIFICATIONS	trocar-flat

ZUSÄTZLICHE INFORMATIONEN

UNIT	Box
TYPE	Kirschner wires
Aesculap	LX014
KLS Martin	22-632-14

AUTHORS

[KIRSCHNER](#)

Medicon

76.15.14

SPECIFICATIONS

[trocar-flat](#)