

KLEPPINGER ELECTRODE; BIPOAR; INSERT

serrated; for Ø 5.0mm shaft

for 628-201-53 and 628-205-53



SKU: 628-248-53

Categories: [Endoscopy](#), [Laparoscopy](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|--------------|---|
| UNIT | pcs. |
| RANGE | HF-Instruments, bipolar |
| TYPE | Electrodes |

ZUSÄTZLICHE INFORMATIONEN

| | |
|--------------|---|
| UNIT | pcs. |
| RANGE | HF-Instruments, bipolar |
| TYPE | Electrodes |