

KOCHER INTESTINAL CLAMP

straight; serrated longitudinal; elastic

25.0cm (10")



SKU: 160-234-25

Categories: [Open Surgery](#), [Stomach](#), [Intestines](#), [Rectum](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Stomach and Intestinal Clamp Forceps
Aesculap	EA151
KLS Martin	30-304-25
AUTHORS	KOCHER
Medicon	50.22.26
SPECIFICATIONS	serrated longitudinal , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Stomach and Intestinal Clamp Forceps
Aesculap	EA151
KLS Martin	30-304-25

AUTHORS

[KOCHER](#)

Medicon

50.22.26

SPECIFICATIONS

[serrated longitudinal](#), [straight](#)