

KOCHER INTESTINAL CLAMP

straight; serrated longitudinal; elastic

28.0cm (11.2")

SKU: 160-234-28

Categories: [Open Surgery](#), [Stomach](#), [Intestines](#), [Rectum](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Stomach and Intestinal Clamp Forceps
Aesculap	EA152
KLS Martin	30-304-28
AUTHORS	KOCHER
Medicon	50.22.28
SPECIFICATIONS	serrated longitudinal , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Stomach and Intestinal Clamp Forceps
Aesculap	EA152
KLS Martin	30-304-28

AUTHORS

[KOCHER](#)

Medicon

50.22.28

SPECIFICATIONS

[serrated longitudinal](#), [straight](#)