

KOCHER INTESTINAL CLAMP

straight; serrated longitudinal; elastic

22.0cm (8.8")

SKU: 160-234-22

Categories: [Open Surgery](#), [Stomach](#), [Intestines](#), [Rectum](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Stomach and Intestinal Clamp Forceps
Aesculap	EA150
KLS Martin	30-304-22
AUTHORS	KOCHER
Medicon	50.22.23
SPECIFICATIONS	serrated longitudinal , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Stomach and Intestinal Clamp Forceps
Aesculap	EA150
KLS Martin	30-304-22

AUTHORS

[KOCHER](#)

Medicon

50.22.23

SPECIFICATIONS

[serrated longitudinal](#), [straight](#)