

## KOCHER INTESTINAL CLAMP

straight; serrated longitudinal; elastic

22.0cm (8.8")



**SKU:** 160-234-22

**Categories:** [Open Surgery](#), [Stomach](#), [Intestines](#), [Rectum](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Stomach and Intestinal Clamp Forceps</a>
<b>Aesculap</b>	EA150
<b>KLS Martin</b>	30-304-22
<b>AUTHORS</b>	<a href="#">KOCHER</a>
<b>Medicon</b>	50.22.23
<b>SPECIFICATIONS</b>	<a href="#">serrated longitudinal</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Stomach and Intestinal Clamp Forceps</a>
<b>Aesculap</b>	EA150
<b>KLS Martin</b>	30-304-22

## AUTHORS

[KOCHER](#)

## Medicon

50.22.23

## SPECIFICATIONS

[serrated longitudinal](#), [straight](#)