


KOCHER-NIPPON HAEMOSTATIC FORCEPS

straight; serrated; 1X2 teeth

16.5cm (6.4")



SKU: 040-406-16
Categories: [Open Surgery](#), [Hemostats](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps , Dissecting and Ligature Forceps
KLS Martin	13-382-16
AUTHORS	KOCHER , NIPPON
SPECIFICATIONS	serrated , straight , toothed



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps , Dissecting and Ligature Forceps
KLS Martin	13-382-16
AUTHORS	KOCHER , NIPPON
SPECIFICATIONS	serrated , straight , toothed