

KOCHER-NIPPON HAEMOSTATIC FORCEPS

straight; serrated; 1X2 teeth

18.5cm (7.2")



SKU: 040-406-18
Categories: [Open Surgery](#), [Hemostats](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps , Dissecting and Ligature Forceps
KLS Martin	13-382-18
AUTHORS	KOCHER , NIPPON
SPECIFICATIONS	serrated , straight , toothed



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps , Dissecting and Ligature Forceps
KLS Martin	13-382-18
AUTHORS	KOCHER , NIPPON
SPECIFICATIONS	serrated , straight , toothed