

KOCHER-OCHSNER HAEMOSTATIC FORCEPS

curved; serrated; 1X2 teeth

30.0cm (12")



SKU: 040-351-30
Categories: [Open Surgery](#), [Hemostats](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps
Aesculap	BH655
KLS Martin	13-415-30
AUTHORS	KOCHER , OCHSNER
SPECIFICATIONS	curved , serrated , toothed



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps
Aesculap	BH655
KLS Martin	13-415-30
AUTHORS	KOCHER , OCHSNER

SPECIFICATIONS

[curved](#), [serrated](#), [toothed](#)