

KOCHER-OCHSNER HAEMOSTATIC FORCEPS

straight; serrated; 1X2 teeth

20.0cm (8")



SKU: 040-350-20
Categories: [Open Surgery](#), [Hemostats](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps
Aesculap	BH646
KLS Martin	13-414-20
AUTHORS	KOCHER , OCHSNER
Medicon	16.12.20
SPECIFICATIONS	serrated , straight , toothed



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps
Aesculap	BH646
KLS Martin	13-414-20

AUTHORS

[KOCHER](#), [OCHSNER](#)

Medicon

16.12.20

SPECIFICATIONS

[serrated](#), [straight](#), [toothed](#)