

KOCHER-OCHSNER HAEMOSTATIC FORCEPS

straight; serrated; 1X2 teeth

30.0cm (12")



SKU: 040-350-30
Categories: [Open Surgery](#), [Hemostats](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps
Aesculap	BH654
KLS Martin	13-414-30
AUTHORS	KOCHER , OCHSNER
SPECIFICATIONS	serrated , straight , toothed



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps
Aesculap	BH654
KLS Martin	13-414-30
AUTHORS	KOCHER , OCHSNER

SPECIFICATIONS

[serrated](#), [straight](#), [toothed](#)