

KÖNIG LIGATURE CONDUCTOR

curved; 5mm


19.5cm (7.8")

SKU: 056-508-05
 Categories: [Open Surgery](#), [Suture](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Ligature Conductors
Aesculap	BM901
KLS Martin	20-875-05
AUTHORS	KOENIG
Medicon	13.32.05
SPECIFICATIONS	curved



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Ligature Conductors
Aesculap	BM901
KLS Martin	20-875-05

AUTHORS

[KOENIG](#)

Medicon

13.32.05

SPECIFICATIONS

[curved](#)