

## KOGAN ENDOSPECULA

without ratchet; 4mm

28.0cm (11")




**SKU:** 180-090-04

**Categories:** [Gynecology](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Endospecula</a>
<b>Aesculap</b>	EM710
<b>KLS Martin</b>	32-381-24
<b>AUTHORS</b>	<a href="#">KOGAN</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Endospecula</a>
<b>Aesculap</b>	EM710
<b>KLS Martin</b>	32-381-24
<b>AUTHORS</b>	<a href="#">KOGAN</a>