

**KOGAN ENDOSPECULA**

with ratchet; 4mm  
24.0cm (9.6"); insulated



**SKU:** 180-093-04  
**Categories:** [Gynecology](#), [Open Surgery](#)

**ZUSÄTZLICHE INFORMATIONEN**

|                |                             |
|----------------|-----------------------------|
| <b>UNIT</b>    | <a href="#">pcs.</a>        |
| <b>TYPE</b>    | <a href="#">Endospecula</a> |
| <b>AUTHORS</b> | <a href="#">KOGAN</a>       |

**ZUSÄTZLICHE INFORMATIONEN**

|                |                             |
|----------------|-----------------------------|
| <b>UNIT</b>    | <a href="#">pcs.</a>        |
| <b>TYPE</b>    | <a href="#">Endospecula</a> |
| <b>AUTHORS</b> | <a href="#">KOGAN</a>       |