

KROENLEIN HAEMOSTATIC FORCEPS

straight; serrated

13.0cm (5.2")



SKU: 040-369-13

Categories: [Open Surgery](#), [Hemostats](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps
Aesculap	BH300
KLS Martin	13-182-13
AUTHORS	KROENLEIN
Medicon	16.06.43
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps
Aesculap	BH300
KLS Martin	13-182-13

AUTHORS

[KROENLEIN](#)

Medicon

16.06.43

SPECIFICATIONS

[serrated](#), [straight](#)