

KRONECKER LIGATURE NEEDLE

right; with round handle

15.0cm (6")

SKU: 056-470-01

Categories: [Open Surgery](#), [Suture](#)



ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|----------------------------------|
| UNIT | pcs. |
| TYPE | Ligature Needles |
| Aesculap | BM801 |
| KLS Martin | 20-835-13 |
| AUTHORS | KRONECKER |
| Medicon | 13.11.14 |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|----------------------------------|
| UNIT | pcs. |
| TYPE | Ligature Needles |
| Aesculap | BM801 |
| KLS Martin | 20-835-13 |
| AUTHORS | KRONECKER |

Medicon

13.11.14