

## LAHEY GALL DUCT CLAMP

curved; serrated

23.0cm (9.2")



**SKU:** 040-501-23

**Categories:** [Hemostats](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Gall Duct Forceps</a>
<b>Aesculap</b>	BJ060
<b>KLS Martin</b>	13-513-23
<b>AUTHORS</b>	<a href="#">LAHEY</a>
<b>Medicon</b>	17.11.23
<b>SPECIFICATIONS</b>	<a href="#">curved</a> , <a href="#">serrated</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Gall Duct Forceps</a>
<b>Aesculap</b>	BJ060
<b>KLS Martin</b>	13-513-23

**AUTHORS**

[LAHEY](#)

**Medicon**

17.11.23

**SPECIFICATIONS**

[curved](#), [serrated](#)