

LANE TOWEL FORCEPS

straight; toothed


15.0cm (6")



SKU: 030-040-15
Categories: [Open Surgery](#), [Hemostats](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Towel Forceps
Aesculap	BF451
AUTHORS	LANE
Medicon	17.59.15
SPECIFICATIONS	straight , toothed



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Towel Forceps
Aesculap	BF451
AUTHORS	LANE
Medicon	17.59.15

SPECIFICATIONS

[straight](#), [toothed](#)