

## LANGE EXCHANGEABLE TIP

fig. 1; oval

for BRÜNINGS interchangeable tubes



**SKU:** 131-488-00

**Categories:** [Open Surgery](#), [Goitre](#), [Laryngology](#), [Mouth](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Universal handles and tips</a>
<b>Aesculap</b>	OP220
<b>KLS Martin</b>	39-624-01
<b>AUTHORS</b>	<a href="#">LANGE</a>
<b>Medicon</b>	73.43.15

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Universal handles and tips</a>
<b>Aesculap</b>	OP220
<b>KLS Martin</b>	39-624-01
<b>AUTHORS</b>	<a href="#">LANGE</a>

**Medicon**

73.43.15