

## LANGE-HOHMANN BONE LEVER

22mm; wide tip with groove

25.0cm (10")


**SKU:** 190-574-22

**Categories:** [Open Surgery](#), [Orthopaedics](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Bone Levers</a>
<b>Aesculap</b>	FK160
<b>KLS Martin</b>	23-549-22
<b>AUTHORS</b>	<a href="#">HOHMANN</a> , <a href="#">LANGE</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Bone Levers</a>
<b>Aesculap</b>	FK160
<b>KLS Martin</b>	23-549-22
<b>AUTHORS</b>	<a href="#">HOHMANN</a> , <a href="#">LANGE</a>