

## LAPAROSCOPE HD+

Ø 5.0mm; 30°; 300 mm

autoclavable

**SKU:** 601-005-30HD

**Categories:** [Endoscopy](#), [Equipment](#), [Laparoscopy](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Telescopes](#)

**SCOPE** [Laparoscopy](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Telescopes](#)

**SCOPE** [Laparoscopy](#)