

## LAPAROSCOPE HD+

Ø 10.0mm; 30°; 425 mm

autoclavable



**SKU:** 601-012-30HD

**Categories:** [Endoscopy](#), [Equipment](#), [Laparoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Telescopes</a>
<b>SCOPE</b>	<a href="#">Laparoscopy</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Telescopes</a>
<b>SCOPE</b>	<a href="#">Laparoscopy</a>