

LAPAROSCOPE HD+2

Ø 10.0mm; 0°; 330 mm

autoclavable

SKU: 601-010-00HD2

Categories: [Endoscopy](#), [Equipment](#), [Laparoscopy](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Telescopes](#)

SCOPE [Laparoscopy](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Telescopes](#)

SCOPE [Laparoscopy](#)