

## LAPAROSCOPE HD+4

Ø 5.5mm; 45°; 300 mm

autoclavable

**SKU:** 601-005-40HD4

**Categories:** [Endoscopy](#), [Equipment](#), [Laparoscopy](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Telescopes](#)

**SCOPE** [Laparoscopy](#)



45°

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Telescopes](#)

**SCOPE** [Laparoscopy](#)