

LAPAROSCOPE HD+4

Ø 10,0mm; 45°; 330 mm

autoclavable

SKU: 601-010-40HD4

Categories: [Endoscopy](#), [Equipment](#), [Laparoscopy](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Telescopes](#)

SCOPE [Laparoscopy](#)



45°

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Telescopes](#)

SCOPE [Laparoscopy](#)