

LASER-WORKING ELEMENT, PASSIVE

For flexible laser-probes up to 0.8 mm

FL 302 mm; for 2.9 mm telescopes/ 22 Fr.

Resectoscopy-systems (630-304-01+630-309-01) Closed handle; LuerLock connector



SKU: 630-046-08

Categories: [Endoscopy](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Instruments for Laser Therapy
TYPE	Working Elements

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Instruments for Laser Therapy
TYPE	Working Elements