

LASER-WORKING ELEMENT

passive cutting; closed handle
 for fibres up to Ø 0.8mm

Compatible with Resectoscopy-Systems 24/26 Fr. (630-204-06 /-26 & 630-209-06/-26) 27/28.5 Fr. (630-204-08/-28 & 630-209-08/-28) Urethrotome sheath 21 Fr. (630-400-21)



SKU: 630-047-08
Categories: [Endoscopy](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Instruments for Laser Therapy
TYPE	Working Elements

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Instruments for Laser Therapy
TYPE	Working Elements